

Product datasheet for **MR211293L4V**

Gaa (NM_001159324) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Gaa (NM_001159324) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Gaa
Synonyms:	E430018M07Rik
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_001159324
ORF Size:	2859 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR211293).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001159324.2
RefSeq Size:	3448 bp
RefSeq ORF:	2862 bp
Locus ID:	14387
UniProt ID:	P70699
Cytogenetics:	11 83.35 cM



[View online »](#)

Gene Summary:

This gene encodes a lysosomal acid glucosidase that is involved in the degradation of glycogen. The encoded preproprotein undergoes proteolytic processing to generate a mature enzyme that cleaves alpha-1-4 and alpha-1-6 glycosidic bonds of glycogen, maltose and intermediate oligosaccharides within the lysosome. Mice lacking the encoded protein exhibit symptoms similar to human Pompe syndrome such as accumulation of glycogen in cardiac and skeletal muscle lysosomes resulting in reduced mobility and strength. Alternative splicing of this gene results in multiple transcript variants. [provided by RefSeq, Nov 2015]