

## Product datasheet for **MR211303L4V**

### Gprasp1 (BC058970) Mouse Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Gprasp1 (BC058970) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Gprasp1
Synonyms:	GASP1, GASP
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	BC058970
ORF Size:	2871 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR211303).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">BC058970</a> , <a href="#">AAH58970</a>
RefSeq Size:	4570 bp
RefSeq ORF:	2873 bp
Locus ID:	67298
Cytogenetics:	X F1
Gene Summary:	Modulates lysosomal sorting and functional down-regulation of a variety of G-protein coupled receptors. Targets receptors for degradation in lysosomes via its interaction with BECN2.[UniProtKB/Swiss-Prot Function]



[View online »](#)