

Product datasheet for **MR211311L3V**

Efcab6 (NM_001161628) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Efcab6 (NM_001161628) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Efcab6
Synonyms:	4931407K02Rik; 4932408N08Rik; bM150J22.2; mKIAA1672
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001161628
ORF Size:	2880 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR211311).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001161628.1
RefSeq Size:	2924 bp
RefSeq ORF:	2487 bp
Locus ID:	77627
UniProt ID:	Q6P1E8
Cytogenetics:	15 E1-E2
Gene Summary:	Negatively regulates the androgen receptor by recruiting histone deacetylase complex, and protein DJ-1 antagonizes this inhibition by abrogation of this complex.[UniProtKB/Swiss-Prot Function]



[View online »](#)