

Product datasheet for MR211315L3

Fgd1 (NM_008001) Mouse Tagged Lenti ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Fgd1 (NM_008001) Mouse Tagged Lenti ORF Clone
Tag:	Myc-DDK
Symbol:	Fgd1
Synonyms:	ZFYVE; ZFYVE3
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
E. coli Selection:	Chloramphenicol (34 ug/mL)
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR211315).
Restriction Sites:	SgfI-MluI
Cloning Scheme:	

Cloning sites used for ORF Shuttling:



* The last codon before the Stop codon of the ORF.

ACCN:	NM_008001
ORF Size:	2880 bp



[View online »](#)

OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_008001.4 , NP_032027.2
RefSeq Size:	3516 bp
RefSeq ORF:	2883 bp
Locus ID:	14163
UniProt ID:	P52734
Cytogenetics:	X 68.46 cM
Gene Summary:	This gene encodes a member of family of Rho-specific guanine nucleotide exchange factors. Rho-specific guanine nucleotide exchange factors catalyze the exchange of GDP for GTP and activate small GTPases, which function as molecular switches in signaling. This protein specifically binds cell division cycle 42, a Rho (Ras homology) GTPase. Investigations in mouse suggest that this protein is important for skeletal mineralization and for regulating the actin cytoskeleton. In humans, mutations in this gene are associated with faciogenital dysplasia, also known as Aarskog-Scott syndrome. [provided by RefSeq, Mar 2014]