

Product datasheet for **MR211329L3V**

Kif5b (NM_008448) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Kif5b (NM_008448) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Kif5b
Synonyms:	AL022807; Khc; Khcs; Kns1; Ukhc
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_008448
ORF Size:	2889 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR211329).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_008448.3 , NP_032474.2
RefSeq Size:	6030 bp
RefSeq ORF:	2892 bp
Locus ID:	16573
UniProt ID:	Q61768
Cytogenetics:	18 4.46 cM



[View online »](#)

Gene Summary:

Microtubule-dependent motor required for normal distribution of mitochondria and lysosomes. May be involved in the mechanisms of growth arrest induced by exposure to DNA-damaging drugs or by cellular senescence (PubMed:9657148). Can induce formation of neurite-like membrane protrusions in non-neuronal cells in a ZFYVE27-dependent manner (PubMed:21976701). Regulates centrosome and nuclear positioning during mitotic entry. During the G2 phase of the cell cycle in a BICD2-dependent manner, antagonizes dynein function and drives the separation of nuclei and centrosomes. Required for anterograde axonal transportation of MAPK8IP3/JIP3 which is essential for MAPK8IP3/JIP3 function in axon elongation (By similarity).[UniProtKB/Swiss-Prot Function]