

## Product datasheet for **MR211332L3V**

### **Rtn3 (NM\_001003934) Mouse Tagged ORF Clone Lentiviral Particle**

#### **Product data:**

Product Type:	Lentiviral Particles
Product Name:	Rtn3 (NM_001003934) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Rtn3
Synonyms:	RTN3-A1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001003934
ORF Size:	2892 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR211332).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_001003934.2</a> , <a href="#">NP_001003934.1</a>
RefSeq Size:	5038 bp
RefSeq ORF:	2895 bp
Locus ID:	20168
UniProt ID:	<a href="#">Q9ES97</a>
Cytogenetics:	19 A



[View online »](#)

**Gene Summary:**

May be involved in membrane trafficking in the early secretory pathway. Inhibits BACE1 activity and amyloid precursor protein processing. May induce caspase-8 cascade and apoptosis. May favor BCL2 translocation to the mitochondria upon endoplasmic reticulum stress (By similarity). Induces the formation of endoplasmic reticulum tubules (PubMed:24262037).[UniProtKB/Swiss-Prot Function]