

## Product datasheet for **MR211345L3V**

### **Ppp1r9a (BC053748) Mouse Tagged ORF Clone Lentiviral Particle**

#### **Product data:**

Product Type:	Lentiviral Particles
Product Name:	Ppp1r9a (BC053748) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Ppp1r9a
Synonyms:	5330407E15, mKIAA1222, A230094E16Rik
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	BC053748
ORF Size:	2907 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR211345).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">BC053748.1</a>
RefSeq Size:	9473 bp
RefSeq ORF:	2909 bp
Locus ID:	243725
Cytogenetics:	6 1.86 cM



[View online »](#)