

## Product datasheet for **MR211361L3V**

### Exoc4 (NM\_009148) Mouse Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Exoc4 (NM_009148) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Exoc4
Synonyms:	C78892; Sec8; Sec8l1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_009148
ORF Size:	2925 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR211361).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_009148.3</a>
RefSeq Size:	3723 bp
RefSeq ORF:	2928 bp
Locus ID:	20336
UniProt ID:	<a href="#">O35382</a>
Cytogenetics:	6 14.59 cM
Gene Summary:	Component of the exocyst complex involved in the docking of exocytic vesicles with fusion sites on the plasma membrane.[UniProtKB/Swiss-Prot Function]



[View online »](#)