

Product datasheet for **MR211427L3V**

Hltf (NM_009210) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Hltf (NM_009210) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Hltf
Synonyms:	AF010600; BC057116; P113; Smarca3; Snf2l3
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_009210
ORF Size:	3009 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR211427).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_009210.2 , NP_033236.2
RefSeq Size:	4955 bp
RefSeq ORF:	3012 bp
Locus ID:	20585
UniProt ID:	Q6PCN7
Cytogenetics:	3 6.15 cM



[View online »](#)

Gene Summary:

Has both helicase and E3 ubiquitin ligase activities. Possesses intrinsic ATP-dependent nucleosome-remodeling activity. This activity may be required for transcriptional activation or repression of specific target promoters (By similarity). These may include the SERPINE1, to which this protein can bind directly. Plays a role in error-free postreplication repair (PRR) of damaged DNA and maintains genomic stability through acting as a ubiquitin ligase for 'Lys-63'-linked polyubiquitination of chromatin-bound PCNA (By similarity).[UniProtKB/Swiss-Prot Function]