

## Product datasheet for **MR211436L3V**

### **Tbc1d9 (BC057027) Mouse Tagged ORF Clone Lentiviral Particle**

#### **Product data:**

Product Type:	Lentiviral Particles
Product Name:	Tbc1d9 (BC057027) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Tbc1d9
Synonyms:	4933431N12Rik; AI847101; AW490653; C76116
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	BC057027
ORF Size:	3018 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR211436).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">BC057027.1</a>
RefSeq Size:	4453 bp
RefSeq ORF:	3020 bp
Locus ID:	71310
Cytogenetics:	8 C2
Gene Summary:	May act as a GTPase-activating protein for Rab family protein(s).[UniProtKB/Swiss-Prot Function]



[View online »](#)