

Product datasheet for **MR211447L4V**

Tcof1 (BC021514) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Tcof1 (BC021514) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Tcof1
Synonyms:	AA408847; AW209012; treacle
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	BC021514
ORF Size:	3039 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR211447).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	BC021514.1
RefSeq Size:	4088 bp
RefSeq ORF:	3041 bp
Locus ID:	21453
Cytogenetics:	18 34.41 cM



[View online »](#)

Gene Summary:

Nucleolar protein that acts as a regulator of RNA polymerase I by connecting RNA polymerase I with enzymes responsible for ribosomal processing and modification. Required for neural crest specification: following monoubiquitination by the BCR(KBTBD8) complex, associates with NOLC1 and acts as a platform to connect RNA polymerase I with enzymes responsible for ribosomal processing and modification, leading to remodel the translational program of differentiating cells in favor of neural crest specification.[UniProtKB/Swiss-Prot Function]