

Product datasheet for **MR211453L3V**

Dennd1a (NM_146122) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Symbol:	Dennd1a
Synonyms:	60304461I9Rik; connecdenn
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_146122
ORF Size:	3048 bp

ORF Nucleotide Sequence: The ORF insert of this clone is exactly the same as(MR211453).

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

RefSeq:	NM_146122.3 , NP_666234.3
RefSeq Size:	4352 bp
RefSeq ORF:	3051 bp
Locus ID:	227801
UniProt ID:	Q8K382
Cytogenetics:	2 B



Gene Summary:

Guanine nucleotide exchange factor (GEF) regulating clathrin-mediated endocytosis through RAB35 activation. Promotes the exchange of GDP to GTP, converting inactive GDP-bound RAB35 into its active GTP-bound form. Regulates clathrin-mediated endocytosis of synaptic vesicles and mediates exit from early endosomes (PubMed:17182770, PubMed:20159556). Binds phosphatidylinositol-phosphates (PtdInsPs), with some preference for PtdIns(3)P (PubMed:20159556).[UniProtKB/Swiss-Prot Function]