

Product datasheet for **MR211476L3V**

Dpyd (NM_170778) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Dpyd (NM_170778) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Dpyd
Synonyms:	AI315208; DPD; E330028L06Rik
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_170778
ORF Size:	3075 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR211476).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_170778.2
RefSeq Size:	4437 bp
RefSeq ORF:	3078 bp
Locus ID:	99586
UniProt ID:	Q8CHR6
Cytogenetics:	3 G1
Gene Summary:	Involved in pyrimidine base degradation. Catalyzes the reduction of uracil and thymine (By similarity).[UniProtKB/Swiss-Prot Function]



[View online »](#)