

## Product datasheet for **MR211522L4V**

### Col16a1 (BC027766) Mouse Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Col16a1 (BC027766) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Col16a1
Synonyms:	2700007F12Rik; A530052M23Rik; AI838662; CA1F
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	BC027766
ORF Size:	3138 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR211522).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">BC027766.1</a>
RefSeq Size:	5160 bp
RefSeq ORF:	3140 bp
Locus ID:	107581
Cytogenetics:	4 D2.2
Gene Summary:	Involved in mediating cell attachment and inducing integrin-mediated cellular reactions, such as cell spreading and alterations in cell morphology.[UniProtKB/Swiss-Prot Function]



[View online »](#)