

Product datasheet for MR211536L3

Trim24 (NM_145076) Mouse Tagged Lenti ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Trim24 (NM_145076) Mouse Tagged Lenti ORF Clone
Tag:	Myc-DDK
Symbol:	Trim24
Synonyms:	A130082H20Rik; AI447469; D430004I05Rik; Tif; TIF1; TIF1-alpha; Tif1a; TIF1alpha
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
E. coli Selection:	Chloramphenicol (34 ug/mL)
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR211536).
Restriction Sites:	SgfI-MluI
Cloning Scheme:	

Cloning sites used for ORF Shuttling:



* The last codon before the Stop codon of the ORF.

ACCN:	NM_145076
ORF Size:	3153 bp



[View online »](#)

OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_145076.4 , NP_659542.3
RefSeq Size:	6191 bp
RefSeq ORF:	3156 bp
Locus ID:	21848
UniProt ID:	Q64127
Cytogenetics:	6 B1
Gene Summary:	The protein encoded by this gene is part of the tripartite-motif containing family (TRIM), which are typified by the RING, B-box type 1, B-box type 2, and coiled-coil region domains. This protein, which also contains a PHD/TTC finger and bromodomain important for regulating nuclear receptors and binding chromatin, has important roles in differentiation, development, and tissue homeostasis. This protein has been reported to regulate the activity of the tumor suppressor p53 and of the retinoic acid receptor. A translocation event between this gene and Braf transforming gene, which results in the fusion protein T18, has been reported in hepatocellular carcinomas. Alternative splicing results in multiple transcript variants that encode different protein isoforms. [provided by RefSeq, Jan 2013]