

## Product datasheet for **MR211552L3V**

### Usp25 (NM\_013918) Mouse Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Usp25 (NM_013918) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Usp25
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_013918
ORF Size:	3165 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR211552).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_013918.2</a> , <a href="#">NP_038946.2</a>
RefSeq Size:	4432 bp
RefSeq ORF:	3168 bp
Locus ID:	30940
UniProt ID:	<a href="#">P57080</a>
Cytogenetics:	16 C3.1



[View online »](#)

**Gene Summary:**

Deubiquitinating enzyme that hydrolyzes ubiquitin moieties conjugated to substrates and thus, functions to process newly synthesized Ubiquitin, to recycle ubiquitin molecules or to edit polyubiquitin chains and prevents proteasomal degradation of substrates. Hydrolyzes both 'Lys-48'- and 'Lys-63'-linked tetraubiquitin chains (By similarity).[UniProtKB/Swiss-Prot Function]