

Product datasheet for **MR211555L4V**

Tcerg1 (BC039185) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Tcerg1 (BC039185) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Tcerg1
Synonyms:	CA150b, p144, ca150, AI428505
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	BC039185
ORF Size:	3171 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR211555).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	BC039185.1
RefSeq Size:	4090 bp
RefSeq ORF:	3173 bp
Locus ID:	56070
Cytogenetics:	18 B3
Gene Summary:	Transcription factor that binds RNA polymerase II and inhibits the elongation of transcripts from target promoters. Regulates transcription elongation in a TATA box-dependent manner (By similarity).[UniProtKB/Swiss-Prot Function]



[View online »](#)