

## Product datasheet for **MR211576L4V**

### Hip1r (NM\_145070) Mouse Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Hip1r (NM_145070) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Hip1r
Synonyms:	AA410023; mKIAA0655
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_145070
ORF Size:	3204 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR211576).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_145070.3</a> , <a href="#">NP_659507.3</a>
RefSeq Size:	4402 bp
RefSeq ORF:	3207 bp
Locus ID:	29816
UniProt ID:	<a href="#">Q9JKY5</a>
Cytogenetics:	5 F



[View online »](#)

**Gene Summary:**

Component of clathrin-coated pits and vesicles, that may link the endocytic machinery to the actin cytoskeleton. Binds 3-phosphoinositides (via ENTH domain). May act through the ENTH domain to promote cell survival by stabilizing receptor tyrosine kinases following ligand-induced endocytosis.[UniProtKB/Swiss-Prot Function]