

Product datasheet for **MR211580L3V**

Ube3b (NM_054093) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Ube3b (NM_054093) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Ube3b
Synonyms:	AI449831; AU020130
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_054093
ORF Size:	3213 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR211580).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_054093.2
RefSeq Size:	5191 bp
RefSeq ORF:	3213 bp
Locus ID:	117146
UniProt ID:	Q9ES34
Cytogenetics:	5 F
Gene Summary:	E3 ubiquitin-protein ligase which accepts ubiquitin from an E2 ubiquitin-conjugating enzyme in the form of a thioester and then directly transfers the ubiquitin to targeted substrates. [UniProtKB/Swiss-Prot Function]



[View online »](#)