

Product datasheet for **MR211631L2V**

Pitpnm2 (BC070452) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Pitpnm2 (BC070452) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Pitpnm2
Synonyms:	NIR3, RDGBA2, mKIAA1457
Mammalian Cell Selection:	None
Vector:	pLenti-C-mGFP (PS100071)
Tag:	mGFP
ACCN:	BC070452
ORF Size:	3294 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR211631).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	BC070452 , AAH70452
RefSeq Size:	6235 bp
RefSeq ORF:	3296 bp
Locus ID:	19679
Cytogenetics:	5 63.39 cM
Gene Summary:	Catalyzes the transfer of phosphatidylinositol and phosphatidylcholine between membranes (in vitro). Binds calcium ions (By similarity).[UniProtKB/Swiss-Prot Function]



[View online »](#)