

## Product datasheet for **MR211657L2V**

### Lrig3 (NM\_177152) Mouse Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Lrig3 (NM_177152) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Lrig3
Synonyms:	9030421L11Rik; 9130004I02Rik; 9430095K15Rik; mKIAA3016
Mammalian Cell Selection:	None
Vector:	pLenti-C-mGFP (PS100071)
Tag:	mGFP
ACCN:	NM_177152
ORF Size:	3351 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR211657).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_177152.5</a> , <a href="#">NP_796126.4</a>
RefSeq Size:	4050 bp
RefSeq ORF:	3354 bp
Locus ID:	320398
UniProt ID:	<a href="#">Q6P1C6</a>
Cytogenetics:	10 D3
Gene Summary:	Plays a role in craniofacial and inner ear morphogenesis during embryonic development. Acts within the otic vesicle epithelium to control formation of the lateral semicircular canal in the inner ear, possibly by restricting the expression of NTN1.[UniProtKB/Swiss-Prot Function]



[View online »](#)