

Product datasheet for **MR211684L4V**

Pcdh18 (NM_130448) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Pcdh18 (NM_130448) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Pcdh18
Synonyms:	3110038E07Rik; BB095589; PCDH68L
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_130448
ORF Size:	3402 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR211684).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_130448.3 , NP_569715.3
RefSeq Size:	5168 bp
RefSeq ORF:	3405 bp
Locus ID:	73173
UniProt ID:	Q8VHR0
Cytogenetics:	3 C
Gene Summary:	Potential calcium-dependent cell-adhesion protein.[UniProtKB/Swiss-Prot Function]



[View online »](#)