

Product datasheet for **MR211755L4V**

Atn1 (NM_007881) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Atn1 (NM_007881) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Atn1
Synonyms:	Atr1; atrophin-1; Drpla
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_007881
ORF Size:	3525 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR211755).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_007881.4 , NP_031907.2
RefSeq Size:	4426 bp
RefSeq ORF:	3528 bp
Locus ID:	13498
UniProt ID:	O35126
Cytogenetics:	6 59.17 cM



[View online »](#)

Gene Summary:

Transcriptional corepressor. Corepressor of MTG8 transcriptional repression. Has some intrinsic repression activity which is independent of the number of the poly-Q repeats (By similarity). Recruits NR2E1 to repress transcription. Promotes vascular smooth cell (VSMC) migration and orientation.[UniProtKB/Swiss-Prot Function]