

Product datasheet for **MR211756L3V**

Atp11b (NM_029570) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Atp11b (NM_029570) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Atp11b
Synonyms:	1110019I14Rik; ATPIF; ATPIR; mKIAA0956
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_029570
ORF Size:	3525 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR211756).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_029570.3
RefSeq Size:	4883 bp
RefSeq ORF:	3528 bp
Locus ID:	76295
UniProt ID:	Q6DFW5
Cytogenetics:	3 B



[View online »](#)

Gene Summary:

Catalytic component of a P4-ATPase flippase complex which catalyzes the hydrolysis of ATP coupled to the transport of aminophospholipids, phosphatidylserines (PS) and phosphatidylethanolamines (PE), from the outer to the inner leaflet of intracellular membranes. May contribute to the maintenance of membrane lipid asymmetry in endosome compartment.[UniProtKB/Swiss-Prot Function]