

Product datasheet for **MR211790L3V**

Sin3a (BC052716) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Sin3a (BC052716) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Sin3a
Synonyms:	AW553200; mKIAA4126; mSin3A; Sin3
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	BC052716
ORF Size:	3591 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR211790).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	BC052716 , AAH52716
RefSeq Size:	4820 bp
RefSeq ORF:	3593 bp
Locus ID:	20466
Cytogenetics:	9 30.89 cM



[View online »](#)

Gene Summary:

Acts as a transcriptional repressor. Corepressor for REST. Interacts with MXI1 to repress MYC responsive genes and antagonize MYC oncogenic activities. Also interacts with MXD1-MAX heterodimers to repress transcription by tethering SIN3A to DNA. Acts cooperatively with OGT to repress transcription in parallel with histone deacetylation. Involved in the control of the circadian rhythms. Required for the transcriptional repression of circadian target genes, such as PER1, mediated by the large PER complex through histone deacetylation. Cooperates with FOXK1 to regulate cell cycle progression probably by repressing cell cycle inhibitor genes expression (PubMed:22476904). Required for cortical neuron differentiation and callosal axon elongation (PubMed:27399968).[UniProtKB/Swiss-Prot Function]