

Product datasheet for **MR211811L4V**

Zeb2 (NM_015753) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Zeb2 (NM_015753) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Zeb2
Synonyms:	9130203F04Rik; D130016B08Rik; SIP1; Zfhx1b; Zfx1b; Zfxh1b
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_015753
ORF Size:	3645 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR211811).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_015753.3
RefSeq Size:	9113 bp
RefSeq ORF:	3648 bp
Locus ID:	24136
UniProt ID:	Q9R0G7
Cytogenetics:	2 B
Gene Summary:	Transcriptional inhibitor that binds to DNA sequence 5'-CACCT-3' in different promoters. Represses transcription of E-cadherin.[UniProtKB/Swiss-Prot Function]



[View online »](#)