

Product datasheet for **MR211821L3V**

Nav1 (BC059840) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Nav1 (BC059840) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Nav1
Synonyms:	9530089B19, POMFIL3, C230080M11Rik
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	BC059840
ORF Size:	3663 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR211821).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	BC059840 , AAH59840
RefSeq Size:	3938 bp
RefSeq ORF:	3665 bp
Locus ID:	215690
Cytogenetics:	1 E4
Gene Summary:	May be involved in neuronal migration.[UniProtKB/Swiss-Prot Function]



[View online »](#)