

## Product datasheet for **MR211834L4V**

### Plcb3 (NM\_008874) Mouse Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Plcb3 (NM_008874) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Plcb3
Synonyms:	mKIAA4098
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_008874
ORF Size:	3702 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR211834).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_008874.4</a> , <a href="#">NP_032900.2</a>
RefSeq Size:	4273 bp
RefSeq ORF:	3705 bp
Locus ID:	18797
UniProt ID:	<a href="#">P51432</a>
Cytogenetics:	19 5.1 cM
Gene Summary:	The production of the second messenger molecules diacylglycerol (DAG) and inositol 1,4,5-trisphosphate (IP3) is mediated by activated phosphatidylinositol-specific phospholipase C enzymes.[UniProtKB/Swiss-Prot Function]



[View online »](#)