

## Product datasheet for **MR211836L4V**

### Phka2 (NM\_001177879) Mouse Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Phka2 (NM_001177879) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Phka2
Synonyms:	6330505C01Rik; D330034O08; Phk
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_001177879
ORF Size:	3705 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR211836).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_001177879.1</a>
RefSeq Size:	4528 bp
RefSeq ORF:	3708 bp
Locus ID:	110094
UniProt ID:	<a href="#">Q8BWJ3</a>
Cytogenetics:	X 73.95 cM
Gene Summary:	Phosphorylase b kinase catalyzes the phosphorylation of serine in certain substrates, including troponin I. The alpha chain may bind calmodulin.[UniProtKB/Swiss-Prot Function]



[View online »](#)