

Product datasheet for **MR211841L3V**

Slc4a3 (BC080271) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Slc4a3 (BC080271) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Slc4a3
Synonyms:	A930038D23Rik; Ae3
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	BC080271
ORF Size:	3717 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR211841).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	BC080271.1
RefSeq Size:	4092 bp
RefSeq ORF:	3719 bp
Locus ID:	20536
Cytogenetics:	1 39.16 cM
Gene Summary:	Plasma membrane anion exchange protein.[UniProtKB/Swiss-Prot Function]


[View online »](#)