

Product datasheet for **MR211851L3V**

Tbc1d9b (NM_029745) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Tbc1d9b (NM_029745) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Tbc1d9b
Synonyms:	2700008N14Rik; AU019384; RP23-319B15.5
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_029745
ORF Size:	3738 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR211851).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_029745.2 , NP_084021.2
RefSeq Size:	5225 bp
RefSeq ORF:	3741 bp
Locus ID:	76795
UniProt ID:	Q5SVR0
Cytogenetics:	11 B1.3
Gene Summary:	May act as a GTPase-activating protein for Rab family protein(s).[UniProtKB/Swiss-Prot Function]



[View online »](#)