

Product datasheet for **MR211854L3V**

TPP2 (BC058239) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Tpp2 (BC058239) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Tpp2
Synonyms:	TPP-2; TppII
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	BC058239
ORF Size:	3747 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR211854).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	BC058239.1
RefSeq Size:	3960 bp
RefSeq ORF:	3749 bp
Locus ID:	22019
Cytogenetics:	1 23.5 cM
Gene Summary:	Component of the proteolytic cascade acting downstream of the 26S proteasome in the ubiquitin-proteasome pathway. May be able to complement the 26S proteasome function to some extent under conditions in which the latter is inhibited (By similarity). Stimulates adipogenesis.[UniProtKB/Swiss-Prot Function]



[View online »](#)