

Product datasheet for MR211898

Nfrkb (NM_172766) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Nfrkb (NM_172766) Mouse Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Nfrkb
Synonyms:	A530090G11Rik
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
ORF Nucleotide Sequence:	>MR211898 representing NM_172766 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGGATTCCTTGGACCACATGCTGACAGATCCCCTTGAACCTGGTCTTGTGGAGATGGTCACAGTACAG
GAATCATGGAGGACTGCCTTCTGGGTGGTACCAGGGTTAGTCTGCCCGAAGACCTTCTGGAGGATCCCGA
GATATTCTTTGATGTCGTCAGCCTGTCCACATGGCAGGAAGTATTAAGTGACTCTCAGCGTGAGCACCTC
CAGCAGTTCCTGCCTCGGTTTCTGCGGACAGTGTGGAGCAGCAGAGGGAGCTCATCTTGCCTATTCA
GTGGGGAGAAGTCCGCTTCGGAACCCGCTGCACATCGCCAGAAAGCTTCCGAGATGGACACTTTAA
CCCTGAGGTGGTCAAGTATCGTCAGCTCTGCTTCAAGTCACAGTACAAGCGGTATCTCAACTCCCAACAA
CAGTATTTCCACCGGCTGCTGAAACAGATTCTTGCTTCCAGGAGTGATCTCTTGGAGATGGCTCGTAGGA
GTGGCCCTGCTCTTCCCTTCCACACAAGCACCATTACCCGTCGCCGAAGCCCCGAGGAGCGGGAGTGGCG
GACCCAGCAGCGTTACTTGAAAGTCTTGAGGGAAGTGAAGGAGGAATGTGGGGACACTGCTCTGTCTCT
GATGAGGAGGATCTCAGCTCATGGCTTCCAAGCTCTCCAGCACGTTCTCCTAGTCTGCGGTGCCCTGA
GGTGGTGCCACGCTTCCACACGGATGAAAAGTGCAGATAAAATAGAAGTGGGGGACAGTGACCT
GAAGTTAATGTTAAAGAAGCACCATGAGAAGCGGAAACATCAACCAGATCACCCAGACCTTTTGACAGGG
GACCTGACCCCTCAGTGACATCATGACTCGAGTAAATGCCGCGAGGAAGGGCTCCCTAGCAGCCTTGATG
ATTTGGCTGTTCTTAAAAAAGGTGAAGGAAAAGGAGGAGAAGAAAAAGAAAAATAAAGTTGATTAA
GTCGGAAGCAGAGGATCTGGCTGAGCCCTAAGCAGTACTGAAGGGGTCCCAACTCTCTCACAGGCCCC
TCTCCCCTGGCAATATCGTCTATCAAGGAGGAACCCCTGGAAGACATCAAGCCTTGCCTTGAATCAATG
AAATATCTCCAGTTTTTCTCTCTTCTGTTAGAGATCTTGCTGCTGGAAGTCAAGGCTAGCCTTCTAT
GCTGGAGGATCGAGTTTTGGACTGGCAGTCTTCTCCAGCCAGCTCCCTCAACAGCTGGTTCTCTGCTGCC
CCCAACTGGGCTGAGTTGGTGTGCTGCTGTCAGTATCTTGCTGGAGAAAGCCGAGCGGTTCTTCTA
GTTTCTCTCCATTTGTTGAATCAAAGAGAAAACCCAGCAGTGGAAATTGCTTGGTCAATCTCAAGATAA
TGAAAAGGAATTAGCTGCTCTTCCACCTGTGGCTGGAACCAAGACCAGGCCTTCTGTAAGGAAAAAT



[View online >](#)

GAAGACAGCTCAGATGCCATGACACCTGTCCCTCGAGTAAGAAGTACTATGTGGTTCGGCCTAGCACAG
GAGAGGAGAAACGGGTTTTTCAAGAACAGGAGCGTTACAGGTATAGCCAACCCCAATAAGGCATTACCTT
TCGCATGCATGGCTTTGAGTCTGTGGTGGGGCCAGTGAAGGGCGTGTGACAAGGAGACCTCCCTCAAC
AAGGCTCGTGAGCATTACTGCTGCGCTCTGACCGGCTGCCTATGTACCATTCTGTCTTTGTTCCGG
ATGCTGCAGCTCGGCTGCCTAATGGAGAAGGCACTCGGGCAGAGATCTGTGAAGTCTCAAGGACTCTCA
GTTTCTGGCTCCAGATGTACCAGCACTCAGGTGAACACAGTGGTGAAGTGGCCACTGGATCGATTGCAT
TATGAGAAAGACCCTTGTGTGAAATATGACATCGGTGGAAGCTGTGGATCTACCTGCATCGTGACCGGA
GTGAGGAAGAGTTTGAACGGATTCAAGCTCAAGCAGCGGCAGCTAAAGCCAGAAAAGCTCTTCAACA
AAAACCAAGCCGCATCCAAGTGAAGTCCAGTAATAAGGAGGGCTCCACGAAGGGCTTAGCGTCTCT
TCCGAGCAGAGCCAGATGAGCCTCAGTACTAGCATGCCGCCACCCAGTTACGCCTGTAACACCCA
CCACGCCAGCATTGCCACCCCATCTCTCCACCCGTGTCAGCGGTGAACAGAAGTGGCTCTAGCAC
GGTCTCTGAGCCAGCCAGTCTAGTTCCAGGTGTGCTTCTGGTGTCTCACCGACAATGCCACAGTAGGG
ACGATGCTTTCCCGCTCCATCCAGACTCCACCCAGTCTCAGGCCACTGCTCGGTTGTGAGCCATT
CCAGCTCAGCAGGACTACCCAGGTTCCGGTGGTAGCCAGCCAGCCTTCTGTGTTTCCAGCAGTC
AGTAGGGCCAGCAGCCACTACCCAGATGCCAGCAGGACCACAGATTCGTGTGCCAGTACTGCTACT
CAAACAAAGTAGTACCCAGGCAGTTATGGCAACAGTCCAGTCAAGGGCAAAGTGCAGCAGCCTCTG
TGCAACGGCTGGACCTGGGCAGACAGGGCTCACAGTGACGAATCTTCCGGCTGCAGTACGCCCCGTGAG
CAAGACAGCCATGAGTTCTCTGGGAACTCTGCTCCAAGTGCCTCCACCAGCCGTTATTAGAACGTC
ACAGGACAGAACATCAAGCAGGTGTGATCACTGGACAGCTTGGTGTGAAGCCCAGACAGGCGCA
GCATTCACACTCACAGCCAAAATTCGGTATCCAGGTAAGGATGTAAGTGCAGCTGCCGCTTCTCCAT
CACCACAGATGCCAAGGGCCAGACGGTTTTGAGAATCACTCCAGACATGATGGCCACATTGGCCAAGTCT
CAGGTTACCACAGTCAAAGTACTCAGGACCTCTTGGAGCAGGAAGTGGCACTGCAGGCAAAGGCATCT
CTGCTACCTTGCATGTTACTTCCAACCTGTCCACGCGGCTGACAGCCCTGCCAAGGCCCTTCCAGCCAG
TGTCCTTCCATCAGCTCCAGCAGGTAAGTACTGTGGTCAAAGTAACTCCTGATCTCAAGCCGACAGAACT
GCAAATTCAGCTTTTTCGCTTGTATGCCAGCTCTGGTGTAAAGTGTGGCAGATCAGAAGGAAAGAACACAG
TGGCCTCTTCAAGCAAAACCTGCCGCCACAATCCGCATCGTGCAGGGGCTGGGAGTGTGCCCCCTAA
AGCAGGCCAGACCATTACTGTTGCAGCACATGCAAGCAAGGAGCCTCTGTTGCTGGAGGGTCTGGAAT
GTCCATTCTTCAACGGTGTCTTACCCAGTATAAATGCCGCTGTGTCAAAGACTGTGGCTGTGGCCTCTG
GGGCAACAAGCACCCCATCAGCATCGGGACTGGAGCCCCACCGTGGCAGAGTTCTGTTAACACTAC
AGTTGTGTCCACATCCCAGTCTGGGAACTGCCTACCAGAATCACAGTTCCTCTTCTGTGATTAGCCAG
CCAATGAAGGGCAAGAGCGTGGTACAGCCCCATCATCAAAGGCAACCTTGGAGCCAATCTCAGTGGC
TGGGTGCAACATCATCCTCACGACCATGCCAGCAGGTACCAAAGTCAATTGCTGGCAATAAGCCAGTGAG
TTTCTCACTGCTCAGCAGTTGACGAGCTTACGCAACAAGTCAAGGCTACACAGGTGCCATCCAGACC
GTCCCGCATCCCATCTGCAACAGGGCACGGCTTCTGGCTCCTCAAAGCAGTGTCCACTGTTGTTGTGA
CCACAGCTCCATCTCTAAACAAGCACCTGAGCAGCAG

ACGCGTACGCGGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >MR211898 representing NM_172766
 Red=Cloning site Green=Tags(s)

MDSL D H M L T D P L E L G P C G D G H S T G I M E D C L L G G T R V S L P E D L L E D P E I F F D V V S L S T W Q E V L S D S Q R E H L
 Q Q F L P R F P A D S V E Q Q R E L I L A L F S G E N F R F G N P L H I A Q K L F R D G H F N P E V V K Y R Q L C F K S Q Y K R Y L N S Q Q
 Q Y F H R L L K Q I L A S R S D L L E M A R R S G P A L P F P H K H H S P S R S P E E R W R T Q Q R Y L K V L R E V K E E C G D T A L S S
 D E E D L S S W L P S S P A R S P S P A V P L R V V P T L S T T D M K T A D K I E L G D S D L K L M L K K H H E K R K H Q P D H P D L L T G
 D L T L S D I M T R V N A G R K G S L A A L Y D L A V L K K K V K E E K K K K I K I K S E A E D L A E P L S S T E G V P T L S Q A P
 S P L A I S S I K E E P L E D I K P C L G I N E I S S S F F S L L L E I L L L E S Q A S L P M L E D R V L D W Q S S P A S S L N S W F S A A
 P N W A E L V L P A L Q Y L A G E S R A V P S S F S P F V E F K E K T Q Q W K L L G S Q D N E K E L A A L F H L W L E T K D Q A F C K E N
 E D S S D A M T P V P R V R T D Y V V R P S T G E E K R V F Q E Q E R Y R Y S Q P H K A F T F R M H G F E S V V G P V K G V F D K E T S L N
 K A R E H S L L R S D R P A Y V T I L S L V R D A A A R L P N G E G T R A E I C E L L K D S Q F L A P D V T S T Q V N T V V S G A L D R L H
 Y E K D P C V K Y D I G R K L W I Y L H R D R S E E F E R I H Q A Q A A A K A R K A L Q Q K P P S K V K S S N K E G S T K G L S G P
 S E Q S Q M S L S D S M P P T P V T P T T P A L P T P I S P P P V S A V N R S G S S T V S E P A Q S S S G V L L V S S P T M P Q L G
 T M L S P A S I Q T P P S S Q A T A R V V S H S S A G L P Q V R V V A Q P S L P A V S Q Q S V G P A Q P L P Q M P A G P Q I R V P V T A T
 Q T K V V P Q A V M A T V P V K G Q T A A A S V Q R P G P G Q T G L T V T N L P A A V S P V S K T A M S S P G N S A P S A S T T A V I Q N V
 T G Q N I I K Q V S I T G Q L G V K P Q T G S S I P L T A T N F R I Q G K D V L R L P P S S I T T D A K G Q T V L R I T P D M M A T L A K S
 Q V T T V K L T Q D L F G A G S T A G K G I S A T L H V T S N P V H A A D S P A K A P S A S V P S S A P A G T T V V K V T P D L K P T E T
 A N S A F R L M P A L G V S V A D Q K G K N T V A S S E A K P A A T I R I V Q L G V M P P K A G Q T I T V A A H A K Q G A S V A G G S G T
 V H S S T V S L P S I N A A V S K T V A V A S G A T S T P I S I G T G A P T V R Q V P V N T T V V S T S Q S G K L P T R I T V P L S V I S Q
 P M K G K S V V T A P I I K G N L G A N L S G L G R N I I L T T M P A G T K L I A G N K P V S F L T A Q Q L Q L Q Q Q G A T Q V R I Q T
 V P A S H L Q Q G T A S G S S K A V S T V V V T T A P S P K Q A P E Q Q

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Chromatograms: https://cdn.origene.com/chromatograms/mm9099_a04.zip

Restriction Sites: SgfI-MluI

Cloning Scheme:




ACCN: NM_172766

ORF Size: 3888 bp

OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_172766.3 , NP_766354.2
RefSeq Size:	5181 bp
RefSeq ORF:	3891 bp
Locus ID:	235134
UniProt ID:	Q6PIJ4
Cytogenetics:	9 A4
MW:	138.8 kDa
Gene Summary:	Binds to the DNA consensus sequence 5'-GGGGAATCTCC-3'. [UniProtKB/Swiss-Prot Function]

Product images:



Circular map for MR211898