

Product datasheet for **MR211901L4V**

Phldb2 (BC050915) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Phldb2 (BC050915) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Phldb2
Synonyms:	LL5b, LL5beta
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	BC050915
ORF Size:	3906 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR211901).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	BC050915.1
RefSeq Size:	5380 bp
RefSeq ORF:	3908 bp
Locus ID:	208177
Cytogenetics:	16 B5
Gene Summary:	Seems to be involved in the assembly of the postsynaptic apparatus. May play a role in acetyl-choline receptor (AChR) aggregation in the postsynaptic membrane.[UniProtKB/Swiss-Prot Function]



[View online »](#)