

Product datasheet for **MR211911L4V**

BC033915 (BC063268) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	BC033915 (BC063268) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	BC033915
Synonyms:	5730525O22Rik; 9030204A07Rik; AI447252; AI846254; AI849445; mKIAA0999; Qsk; SIK-3
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	BC063268
ORF Size:	3933 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR211911).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	BC063268.1
RefSeq Size:	6099 bp
RefSeq ORF:	3935 bp
Locus ID:	70661
Cytogenetics:	9 A5.2



[View online »](#)

Gene Summary:

Positive regulator of mTOR signaling that functions by triggering the degradation of DEPTOR, an mTOR inhibitor (By similarity). Required for chondrocyte hypertrophy during skeletogenesis (PubMed:22318228). Negatively regulates cAMP signaling pathway possibly by acting on CRTC2/TORC2 and CRTC3/TORC3 (By similarity). Prevents HDAC4 translocation to the nucleus (PubMed:22318228).[UniProtKB/Swiss-Prot Function]