

## Product datasheet for **MR211928L3V**

### **Kdm3a (NM\_001038695) Mouse Tagged ORF Clone Lentiviral Particle**

#### **Product data:**

|                           |  |
|---------------------------|--|
| Product Type:             | Lentiviral Particles   |
| Product Name:             | Kdm3a (NM_001038695) Mouse Tagged ORF Clone Lentiviral Particle  |
| Symbol:                   | Kdm3a  |
| Synonyms:                 | 1700105C21Rik; C230043E16Rik; JHDM2a; Jmjd1; Jmjd1a; KDM2A; TGSA; Tsga   |
| Mammalian Cell Selection: | Puromycin  |
| Vector:                   | pLenti-C-Myc-DDK-P2A-Puro (PS100092)   |
| Tag:                      | Myc-DDK  |
| ACCN:                     | NM_001038695   |
| ORF Size:                 | 3969 bp  |
| ORF Nucleotide Sequence:  | The ORF insert of this clone is exactly the same as(MR211928).   |
| OTI Disclaimer:           | The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a> |
| OTI Annotation:           | This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.   |
| RefSeq:                   | <a href="#">NM_001038695.2</a> , <a href="#">NP_001033784.2</a>  |
| RefSeq Size:              | 4829 bp  |
| RefSeq ORF:               | 3972 bp  |
| Locus ID:                 | 104263   |
| UniProt ID:               | <a href="#">Q6PCM1</a>   |
| Cytogenetics:             | 6 C1   |



[View online »](#)

**Gene Summary:**

Histone demethylase that specifically demethylates 'Lys-9' of histone H3, thereby playing a central role in histone code. Preferentially demethylates mono- and dimethylated H3 'Lys-9' residue, with a preference for dimethylated residue, while it has weak or no activity on trimethylated H3 'Lys-9'. Demethylation of Lys residue generates formaldehyde and succinate. Involved in hormone-dependent transcriptional activation, by participating in recruitment to androgen-receptor target genes, resulting in H3 'Lys-9' demethylation and transcriptional activation (By similarity). Involved in spermatogenesis by regulating expression of target genes such as PRM1 and TNP1 which are required for packaging and condensation of sperm chromatin (PubMed:17943087). Involved in obesity resistance through regulation of metabolic genes such as PPARA and UCP1.[UniProtKB/Swiss-Prot Function]