

Product datasheet for **MR211929L4V**

Ikbkap (NM_026079) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Ikbkap (NM_026079) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Ikbkap
Synonyms:	3110040G09Rik; 6030413P05; C78473; Elp1; IKAP
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_026079
ORF Size:	4002 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR211929).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_026079.1
RefSeq Size:	6160 bp
RefSeq ORF:	4002 bp
Locus ID:	230233
UniProt ID:	Q7TT37
Cytogenetics:	4 31.66 cM
Gene Summary:	May act as a scaffold protein that may assemble active IKK-MAP3K14 complexes (IKKA, IKKB and MAP3K14/NIK).[UniProtKB/Swiss-Prot Function]



[View online »](#)