

Product datasheet for **MR211933L4V**

Dhx57 (NM_198942) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Symbol:	Dhx57
Synonyms:	AW494914
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_198942
ORF Size:	4005 bp

ORF Nucleotide Sequence: The ORF insert of this clone is exactly the same as(MR211933).

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

RefSeq:	NM_198942.2 , NP_945180.2
RefSeq Size:	4935 bp
RefSeq ORF:	4008 bp
Locus ID:	106794
UniProt ID:	Q6P5D3
Cytogenetics:	17 E3



Gene Summary: Probable ATP-binding RNA helicase.[UniProtKB/Swiss-Prot Function]