

Product datasheet for MR211935L4

Erbb3 (BC106091) Mouse Tagged Lenti ORF Clone

Product data:

Product Type: Expression Plasmids

Product Name: Erbb3 (BC106091) Mouse Tagged Lenti ORF Clone

Tag: mGFP Symbol: Erbb3

Synonyms: C76256; Erbb-3; Erbb3r; Her3

Mammalian Cell Puromycin

Selection:

Vector: pLenti-C-mGFP-P2A-Puro (PS100093)

E. coli Selection: Chloramphenicol (34 ug/mL)


ORF Nucleotide The ORF insert of this clone is exactly the same as(MR211935).

Sequence:

Restriction Sites: Sgfl-Mlul

Cloning Scheme:





^{*} The last codon before the Stop codon of the ORF.




OriGene Technologies, Inc. 9620 Medical Center Drive, Ste 200

CN: techsupport@origene.cn

Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com



Plasmid Map:



ACCN: BC106091 **ORF Size:** 4017 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of

reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing

variants is recommended prior to use. More info

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression

varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube

containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method: 1. Centrifuge at 5,000xg for 5min.

2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.

3. Close the tube and incubate for 10 minutes at room temperature.

4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid

at the bottom.

5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of

shipping when stored at -20°C.

RefSeq: <u>BC106091.1</u>

RefSeq Size: 4911 bp

RefSeq ORF: 4019 bp



Erbb3 (BC106091) Mouse Tagged Lenti ORF Clone - MR211935L4

Locus ID: 13867

Cytogenetics: 10 77.1 cM

Gene Summary: Tyrosine-protein kinase that plays an essential role as cell surface receptor for neuregulins.

Binds to neuregulin-1 (NRG1) and is activated by it; ligand-binding increases phosphorylation

on tyrosine residues and promotes its association with the p85 subunit of

phosphatidylinositol 3-kinase. May also be activated by CSPG5. Involved in the regulation of

myeloid cell differentiation.[UniProtKB/Swiss-Prot Function]