

Product datasheet for **MR211960L4V**

Nup155 (NM_133227) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Nup155 (NM_133227) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Nup155
Synonyms:	D930027M19Rik; mKIAA0791
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_133227
ORF Size:	4173 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR211960).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_133227.3 , NP_573490.3
RefSeq Size:	6420 bp
RefSeq ORF:	4176 bp
Locus ID:	170762
UniProt ID:	Q99P88
Cytogenetics:	15 A1
Gene Summary:	Essential component of nuclear pore complex. Could be essential for embryogenesis. Nucleoporins may be involved both in binding and translocating proteins during nucleocytoplasmic transport.[UniProtKB/Swiss-Prot Function]



[View online »](#)