

Product datasheet for **MR211995L2V**

Bcl9l (BC082304) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Bcl9l (BC082304) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Bcl9l
Synonyms:	MGC7014, DLNB11, BCL9-2, B9L
Mammalian Cell Selection:	None
Vector:	pLenti-C-mGFP (PS100071)
Tag:	mGFP
ACCN:	BC082304
ORF Size:	4371 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR211995).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	BC082304 , AAH82304
RefSeq Size:	5472 bp
RefSeq ORF:	4373 bp
Locus ID:	80288
Cytogenetics:	9 A5.2
Gene Summary:	Transcriptional regulator that acts as an activator. Promotes beta-catenin transcriptional activity. Plays a role in tumorigenesis. Enhances the neoplastic transforming activity of CTNNB1.[UniProtKB/Swiss-Prot Function]



[View online »](#)