

Product datasheet for **MR212008L4V**

Col2a1 (NM_031163) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Col2a1 (NM_031163) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Col2a1
Synonyms:	Co; Col2; Col2a; Col2a-1; Del; Del1; Dmm; L; Lpk; M100413; Rgsc4; Rgsc8; Rgsc413; Rgsc856
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_031163
ORF Size:	4464 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR212008).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_031163.3 , NP_112440.2
RefSeq Size:	5149 bp
RefSeq ORF:	4464 bp
Locus ID:	12824
UniProt ID:	P28481
Cytogenetics:	15 53.97 cM



[View online »](#)

Gene Summary:

This gene encodes the alpha-1 subunit of the fibril-forming type II collagen, the major component of cartilage and the vitreous humor of the eye. The encoded preproprotein forms homotrimeric, triple helical procollagen that undergoes proteolytic processing during fibril formation. Mice harboring certain mutations in this gene exhibit severe chondrodysplasia characterized by short limbs and trunk, craniofacial deformities and cleft palate. A complete lack of the encoded protein in mice results in postnatal lethality. Alternative splicing results in multiple transcript variants encoding different isoforms that may undergo similar proteolytic processing. [provided by RefSeq, Dec 2015]