

## Product datasheet for **MR212015L3V**

### Jarid1c (BC043096) Mouse Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Jarid1c (BC043096) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Jarid1c
Synonyms:	D930009K15Rik; Jarid1c; mKIAA0234; Smcx
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	BC043096
ORF Size:	4530 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR212015).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">BC043096.1</a>
RefSeq Size:	5361 bp
RefSeq ORF:	4532 bp
Locus ID:	20591
Cytogenetics:	X 68.46 cM



[View online »](#)

**Gene Summary:**

Histone demethylase that specifically demethylates 'Lys-4' of histone H3, thereby playing a central role in histone code. Does not demethylate histone H3 'Lys-9', H3 'Lys-27', H3 'Lys-36', H3 'Lys-79' or H4 'Lys-20'. Demethylates trimethylated and dimethylated but not monomethylated H3 'Lys-4'. Participates in transcriptional repression of neuronal genes by recruiting histone deacetylases and REST at neuron-restrictive silencer elements (By similarity). Represses the CLOCK-ARNTL/BMAL1 heterodimer-mediated transcriptional activation of the core clock component PER2.[UniProtKB/Swiss-Prot Function]