

Product datasheet for **MR212020L4V**

Nup214 (BC039282) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Nup214 (BC039282) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Nup214
Synonyms:	CAN, RP23-55I14.4
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	BC039282
ORF Size:	4587 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR212020).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	BC039282.1
RefSeq Size:	4962 bp
RefSeq ORF:	4589 bp
Locus ID:	227720
Cytogenetics:	2 21.97 cM
Gene Summary:	May serve as a docking site in the receptor-mediated import of substrates across the nuclear pore complex.[UniProtKB/Swiss-Prot Function]



[View online »](#)