

Product datasheet for **MR212061L4V**

Pprc1 (BC066048) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Pprc1 (BC066048) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Pprc1
Synonyms:	MGC12121, PRC, mKIAA0595, A930018M19
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	BC066048
ORF Size:	4932 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR212061).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	BC066048.1
RefSeq Size:	5117 bp
RefSeq ORF:	4934 bp
Locus ID:	226169
Cytogenetics:	19 C3
Gene Summary:	Acts as a coactivator during transcriptional activation of nuclear genes related to mitochondrial biogenesis and cell growth. Involved in the transcription coactivation of CREB and NRF1 target genes (By similarity).[UniProtKB/Swiss-Prot Function]



[View online »](#)