

Product datasheet for MR212197

Olfr38 (NM_146986) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Olfr38 (NM_146986) Mouse Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Olfr38
Synonyms:	18A; MOR257-2
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
ORF Nucleotide Sequence:	>MR212197 representing NM_146986 Red=Cloning site Blue=ORF Green=Tags(s)

TTTGTGAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCCGCGATCGCC

ATGGAACCAACAATGAGACACAGCTGCATGAATTCATCCTTCTTGGCCTTTCAGTGACTGGGACTCTC
 AAGTCTCCCTCTTCGTCCTGTTCTTGCTGATGTACCTGGTGACAGTATTGGGAACTTCCTCATATTGT
 TCTGATCAGACTGGACAGCAGACTCCACACTCCCATGTATTTCTTCTCACCAACCTGTCCTGGTGGAT
 GCATCTTATGCCACAAGCATAGTCCCCAGCTGCTGGCTCATTTCTTGGCACACAAAAGCAATTCCAT
 TTCTGAGCTGTGCAGCCAGTTATTTTTTCCCTGGGCTTGGGTGGGATTGAGTTCCTTCTGCTGGCAGT
 GATGGCCTATGACCGCTATGTGGCAGTGTGTGATCCTCTGAGGTATTCAGTCATCATGCATGCAGGGCTC
 TGCACAAGGCTGGTCATCACATCATGGTCAGTGGCTCCATCAACTCTCTTGTGCACACTGCCATCACCT
 TCCAGTGGCCATGTGCACAAATAAATATATTGATCACATAGCCTGTGAAACCCTTGCTGTCGTCAGATT
 GGCTTGTGTGGATACCTCATCTAACAAGATTGCAATCATGGTTTCTAGCATTGTCTTGCTTATGACACCT
 TTCTTCTAGTTCTTCTGTCCTACATCCAGATCATCTCCACCATCTTAAAGATACAGTCCACAGAGGGAA
 GGTGGAAGCCTTTCACACTTGTGCCTCTCACCTTACTATGGTTGTCTGTGCTATGGCATGGCCATTTT
 CACATACATCAACCCCACTCCAGCCCCCTGTCTTCTCAGGAGAAGTTGATGTCTCTGTTCTATGCTATC
 TTGACACCCATGCTGAACCCCATGATTTACAGCCTAAGGAATAAAGAGGTTAAGGGGGCATGGCAGAAAT
 TATTAGGAAATTCTCTGGGTTTGCATCAAAATTGGCCACT

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA


[View online »](#)

Protein Sequence: >MR212197 representing NM_146986
 Red=Cloning site Green=Tags(s)

METNNETQLHEFILLGLSSDWDSQVSLFVLFLMYLVTVLGNFLIIVLIRLDSRLHTPMYFFLTNLSLVD
 ASYATSIVPQLLAHFLATHKAIPFLSCAAQLFFSLGLGGIEFLLAVMAYDRYVAVCDPLRYSVIMHAGL
 CTRLVITSWVSGSINSLVHTAITFQLPMCTNKYIDHIACETLAVVRLACVDTSSNKIAIMVSSIVLLMTP
 FFLVLLSYIQIISTILKIQSTEGRWKAFHTCASHLTMVVL CYGMAIFTYIQPHSSPSVLQEKLMSLFYAI
 LTPMLNPMIYSLRNKEVKGAWQKLLGKFSGFASKLAT



TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites:

SgfI-MluI

Cloning Scheme:

Cloning sites used for ORF Shuttling:



* The last codon before the Stop codon of the ORF

ACCN: NM_146986

ORF Size: 951 bp


OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:	1. Centrifuge at 5,000xg for 5min.
	2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
	3. Close the tube and incubate for 10 minutes at room temperature.
	4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
	5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
Note:	Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.
RefSeq:	<u>NM_146986.1, NP_667197.1</u>
RefSeq Size:	954 bp
RefSeq ORF:	954 bp
Locus ID:	258988
Cytogenetics:	6 20.89 cM
MW:	35.9 kDa
Gene Summary:	Olfactory receptors interact with odorant molecules in the nose, to initiate a neuronal response that triggers the perception of a smell. The olfactory receptor proteins are members of a large family of G-protein-coupled receptors (GPCR) arising from single coding-exon genes. Olfactory receptors share a 7-transmembrane domain structure with many neurotransmitter and hormone receptors and are responsible for the recognition and G protein-mediated transduction of odorant signals. The olfactory receptor gene family is the largest in the genome. The nomenclature assigned to the olfactory receptor genes and proteins for this organism is independent of other organisms. [provided by RefSeq, Jul 2008]

Product images:



Circular map for MR212197