

Product datasheet for MR212374

Mup6 (NM_001081285) Mouse Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: Mup6 (NM_001081285) Mouse Tagged ORF Clone
Tag: Myc-DDK
Symbol: Mup6
Synonyms: 620807; Gm12544; Mup2; OTTMUSG00000007423
Mammalian Cell Selection: Neomycin
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
ORF Nucleotide Sequence: >MR212374 representing NM_001081285
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGAAGCTGCTGCTGTGTTGGGACTGATCCTAGTCTGTGTCCATGCAGAAGAAGCTAGTTCTATGG
GAAGGAACCTTAAATGTAGAAAAGATTAATGGGGAGTGGTATACTATTATCCTGGCCTCTGACAAAAGAGC
AAAGATAGAAGAACATGGCATCATGAGACTTTTTGTGGAGCACATTCATGTCTGGAGAATTCCTTAGGT
TTTAAATCCACTGTTATAGATGAAGAGTGTCCGAAATATTTTTGGTTGCTGACAAAACAGAAAAGG
CTGGTGAATATTCTGTGACATATGATGGATTCAAAAAATTTACAGTACTTAAGACAGACTATGATAATTA
TATTATGTTTCATCTCATTAAATGAAATGAATGGGGAAACCTCCAGCTGATGTCGCTCTATGGCCGAGAA
CCAGATTTGAATTCAGATATCAAAGAAAAGTTTGTAAAATTTGTGAGGAGCATGGAATCATTAGAGAGA
ATATCATTGACTTCACCAAAACCAATCGCTGCCTCCAGGCCCGAGAA

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >MR212374 representing NM_001081285
Red=Cloning site Green=Tags(s)

MKLLLLLCLGLILVCVHAEAEASSMGRNFNVEKINGEYWTIILASDKRAKIEEHGIMRLFVEHIHVLNSLG
FKFHTVIDEECEIFLVADKTEKAGEYSVTYDGFKKFTVLKTDYDNYIMFHLINEMNGETFQLMSLYGRE
PDLNSDIKEKFKVLC EEHGI IRENIIDF TKTNRCLQARE

TRTRPLEQKLI SEEDLAANDILDYKDDDDKV

Chromatograms: https://cdn.origene.com/chromatograms/mm9055_h01.zip



Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_001081285

ORF Size: 537 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001081285.2](#)

RefSeq Size: 540 bp


RefSeq ORF: 540 bp

Locus ID: 620807

Cytogenetics: 4

MW: 20.8 kDa

Product images:



Circular map for MR212374