

Product datasheet for **MR212448**

Szrd1 (NM_001025608) Mouse Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: Szrd1 (NM_001025608) Mouse Tagged ORF Clone
Tag: Myc-DDK
Symbol: Szrd1
Synonyms: 1110022I03; D4Ertd22e
Mammalian Cell Selection: Neomycin
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
ORF Nucleotide Sequence: >MR212448 representing NM_001025608
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGGATCGCC**

ATGGAAGATGAGGAGGTCGCTGAGAGCTGGGAGGAGCGGCAGACAGCGGGAAATAGACAGACGGTTGG
AAAAAACTGAAGATCACACAAAAGAGAGCAGGAAATCCAAATCTCTCCCAAAGTGCCATTGTGAT
TCAAGACGATAGCCTTCCCACGGGGCCCTCCACAGATCCGCATCTCAAGAGGCCACCAGCAACGGT
GTGGTCAGCAGCCCACTCCACCAGCAGGCCAGCCCTTCTGTCAAGTCCCTAGCACAGCGGGAGGCAG
AGTATGCAGAGGCTCGGAGACGGATCCTAGGCAGTGCCAGCCCTGAAGAGGAGCAGGAGAAACCCATCCT
GGACAGGCCAACCAGGATCTCCAACCCGAAGACAGCAGGCAGCCAGCAATGTGATCAGACAGCCGTTG
GGTCTGACGGGTACAAGGCTTCAAACAGCGCAGA

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >MR212448 representing NM_001025608
Red=Cloning site Green=Tags(s)

MEDEEVAESWEEAADSGEIDRRLEKCLKITQKESRKSPPKVPVIVIQDDSLPTGPPPQIRILKRPTSNG
VYSSPNSTRPALPVKSLAQREAEYAEARRRILGSASPEEEQEKPIIDRPTRISQPEDSRQSPNVIRQPL
GPDGSQGFQRR

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Chromatograms: https://cdn.origene.com/chromatograms/mm9064_b12.zip



[View online >](#)

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_001025608

ORF Size: 456 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001025608.2](#), [NP_001020779.1](#)


RefSeq Size: 3297 bp

RefSeq ORF: 459 bp

Locus ID: 213491

UniProt ID: Q6NXN1
 Cytogenetics: 4 73.38 cM
 MW: 17 kDa

Product images:



Circular map for MR212448