

## Product datasheet for **MR212502L4V**

### **Ncoa7 (NM\_172495) Mouse Tagged ORF Clone Lentiviral Particle**

#### **Product data:**

Product Type:	Lentiviral Particles
Product Name:	Ncoa7 (NM_172495) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Ncoa7
Synonyms:	9030406N13Rik
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_172495
ORF Size:	2829 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR212502).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_172495.5</a>
RefSeq Size:	5134 bp
RefSeq ORF:	2832 bp
Locus ID:	211329
UniProt ID:	<a href="#">Q6DFV7</a>
Cytogenetics:	10 A4
Gene Summary:	Enhances the transcriptional activities of several nuclear receptors. Involved in the coactivation of different nuclear receptors, such as ESR1, THRB, PPARG and RARA (By similarity).[UniProtKB/Swiss-Prot Function]



[View online »](#)