

Product datasheet for MR212523

Olfr129 (NM_146327) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Olfr129 (NM_146327) Mouse Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Olfr129
Synonyms:	MOR263-9
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
ORF Nucleotide Sequence:	>MR212523 representing NM_146327 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCCGCGATCGCC

ATGAGTGTCAACTGCTCTCTGTGGCAGGAGAACAGCTTGTCTGTCAAACGCTTTGCATTGCCAAGTTCT
 CTGAGGTCCCTGGAGAATGCTTCCTCCTGTTTACCCTCATCCTCCTCATGTTCTTAGTATCACTGACAGG
 AAATGCACTCATAACCCTTGCCATCTGCACCAGTCCAGCCCTACACACTCCCATGTAATTCTTTCTGGCC
 AACTTGTCTCTCCTGGAGATTGGCTACACTTGTACTGTCATACCGAAGATGTTGCAGAGCCTTGTGAGTG
 AGGCCCGAGGGATCTCTCGGGAGGGATGTGCCACACAGATGTTTTTCTTCACTATTTGGTATAACTGA
 GTGCTGCCTACTGGCAGCCATGGCCTTTGACCGCTGCATGGCCATATGCTCCCCACTCCACTATACAACC
 CGAATGAATCGTGGGGTATGTGCCCATTTGGCAATAGTTTCATGGGGAACAGGTTGCATTGTAGGGTTGG
 GACAAACCAATTTTATTTTCTCCTTAACTTCTGTGGCCCTGTGAGATAGACCACTTCTTCTGTGACCT
 TCCACCTGTCTGGCACTTGCTGTGGAGATACATCCCAAAATGAGGCTGCCATCTTTGTGGTAGCCATC
 CTCTGTATTTCTAGTCCATTTTACTCATCTTATTCCTATGTCCGAATCTGGTTGCAGTGCTGGTGA
 TGCCTTCAACCGAGGGGCGCCATAAAGCTTCTCTACCTGTTCTTCCATCTACTGTAGTCACACTGTT
 TTTTGGCTCAGGATCTATTAACCTACTTGAGGCCAAGTCTAGCCACTCACCAGGAATGGACAACTCTTG
 GCCCTTTTCTACACAGCGGTGACTTCCATGCTGAACCCCATCATCTATAGCTTAAGGAACAAGGAAGTGA
 AGGCAGCACTGAGAAGAACTCTGGGCCTGAAGAAAGTTCTGACTATGAATAGG

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA


[View online »](#)

Protein Sequence: >MR212523 representing NM_146327

Red=Cloning site Green=Tags(s)

MSVNCSLWQENSLSVKRF AFAKFSEVPGECFLLFTLLMFLVSLTGNALITLAICTSPALHTPMYFFLA
NL SLEIGYTCVTPKMLQSLVSEARGISREGCATQMFFFTLFGITECCLAAMAFDRCMAICSPLHYTT
RMNRGVCAHLAIVSWGTCIVGLGQTNFIFSLNFCGCPCEIDHFFCDLPPVLALACGDTSQNEAAIFVVAI
LCISSPFLLLIISYVRILVAVLVMPSPEGRHKAFSTCSSHLLVTLTFFGSGSINYLRPKSSHSHPGMDKLL
ALFYTAVTSMLNP I IYSRLNKEVPAALRRLTLGLKKVLTMMNR



TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Chromatograms: https://cdn.origene.com/chromatograms/mm9083_f01.zip

Restriction Sites: Sgfl-MluI

Cloning Scheme:

Cloning sites used for ORF Shuttling:



* The last codon before the Stop codon of the ORF

ACCN: NM_146327

ORF Size: 963 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

Note: Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.

RefSeq: [NM_146327.2](#), [NP_666439.2](#)

RefSeq Size: 1238 bp

RefSeq ORF: 966 bp


Locus ID: 258324

Cytogenetics: 17 B1

MW: 35.3 kDa

Gene Summary: Olfactory receptors interact with odorant molecules in the nose, to initiate a neuronal response that triggers the perception of a smell. The olfactory receptor proteins are members of a large family of G-protein-coupled receptors (GPCR) arising from single coding-exon genes. Olfactory receptors share a 7-transmembrane domain structure with many neurotransmitter and hormone receptors and are responsible for the recognition and G protein-mediated transduction of odorant signals. The olfactory receptor gene family is the largest in the genome. The nomenclature assigned to the olfactory receptor genes and proteins for this organism is independent of other organisms. [provided by RefSeq, Jul 2008]

Product images:



Circular map for MR212523